**VIASURE**

Real Time PCR Detection Kit

**SARS-CoV-2 & UK Variant  
(S UK, ORF1ab and N genes)**

---

CE IVD

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## Content

1.	Intended use.....	5
2.	Summary and Explanation .....	5
3.	Principle of the procedure .....	6
4.	Reagents provided .....	7
5.	Reagents and equipment to be supplied by the user .....	7
6.	Transport and storage conditions.....	8
7.	Precautions for users .....	8
8.	Test procedure.....	9
9.	Result interpretation .....	10
10.	Limitations of the test .....	13
11.	Quality control .....	14
12.	Performance characteristics.....	15
	ANNEX 1 .....	18
A1.1	Principle of the procedure .....	18
A1.2	Reagents provided.....	18
A1.3	Test procedure .....	19
	ANNEX 2 .....	21
A2.1	Principle of the procedure .....	21
A2.2	Reagents provided.....	21
A2.3	Test procedure .....	22

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## **ENGLISH**

### **1. Intended use**

VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit is a real-time RT-PCR test designed for the qualitative detection of RNA from the SARS-CoV-2 and the HV 69/70 deletion of the *S* gene for SARS-CoV-2 associated to the SARS-CoV-2 VOC202012/01 (lineage B.1.1.7) variant and other variants in nasopharyngeal swabs from individuals suspected of COVID-19 infection by their healthcare professional (HCP). This test is intended for use as an aid in the diagnosis of SARS-CoV-2 as well as variants that carry the HV 69/70 deletion in combination with clinical and epidemiological risk factors. RNA is extracted from respiratory specimens, amplified using RT-PCR and detected using fluorescent reporter dye probes specific for SARS-CoV-2 and HV 69/70 deletion.

### **2. Summary and Explanation**

Coronavirus are enveloped non-segmented positive-sense RNA viruses and belong to Coronaviridae family. There are six coronavirus species known to cause human diseases. Four viruses (229E, OC43, NL63 and HKU1) cause common cold symptoms and the other two (severe acute respiratory syndrome coronavirus (SARS-CoV) and Middle East respiratory syndrome coronavirus (MERS-CoV)) are zoonotic and producing more severe complications. SARS-CoV and MERS-CoV have caused more than 10,000 cumulative cases in the past two decades, with mortality rates of 34% MERS-CoV and 10% SARS-CoV.

In December 2019, some people that worked at or lived around the Huanan seafood market in Wuhan, Hubei Province, China, have presented pneumonia of unknown cause. Deep sequencing analysis of the respiratory samples indicated a novel coronavirus, which was named firstly 2019 novel coronavirus (2019-nCoV) and lately SARS-CoV-2.

Human-to-human transmission of the SARS-CoV-2 has been confirmed, even in the incubation period without symptoms, and the virus causes severe respiratory illness like those SARS-CoV produced. Although the pneumonia is the principal illness associated, a few patients have developed severe pneumonia, pulmonary edema, acute respiratory distress syndrome, or multiple organ failure and death. Center of Disease Control and Prevention (CDC) believes that symptoms of SARS-CoV-2 may appear in as few as 2 days or as long as 14 days after exposure, being the most common fever, cough, myalgia and dyspnea. Less common symptoms are sore throat, headache, diarrhea and vomiting. It seems that people above 60 years old, males, and people with comorbidities most often have severe disease.

Diagnosis of SARS-CoV-2 is performed detecting conventional causes of pneumonia early and detected by next-generation sequencing or real-time RT-PCR methods. Several assays that detect the SARS-CoV-2 have been available, such as China CDC (gene targets, *ORF1ab* and *N*), Charité – Germany (gene targets, *RdRP* and *E*) or US CDC (two targets in *N* gene).

CDC recommends upper respiratory tract specimens (nasopharyngeal (NP) swab, oropharyngeal (OP) swabs, nasal mid-turbinate swab, nasal swab, nasopharyngeal wash/aspirate or nasal wash/aspirate (NW) and saliva specimens collected mainly by a healthcare professional) and/or lower respiratory specimens (sputum, endotracheal aspirate, or bronchoalveolar lavage in patients with more severe respiratory disease) for the

---

identification of SARS-CoV-2. In addition, other clinical specimens as blood, urine and stool may be collected to monitor the presence of the virus.

Since the initial genomic characterization of SARS-CoV-2, the virus has been divided into different genetic groups or clusters (S, L, V, G with GH and GR subgroups). The appearance of mutations is a natural and expected event within the evolution process of the virus. In fact, some specific mutations define the viral genetic groups that are currently circulating globally. The mutations identified to date remain within the expected patterns for a coronavirus. Viruses classified in genetic group G are the most frequent worldwide. Thanks to the genetic sequencing of the pathogen worldwide, it has been possible to establish patterns of dispersal and evolution of the virus.

On December 14, 2020, the United Kingdom declared an increase in the incidence of SARS-CoV-2 in some regions of its country associated with a new variant of the virus with a supposed greater transmission capacity. This variant, called VOC202012/01 (B.1.1.7) presented 23 different mutations: 13 non-synonymous, including a series of mutations in the spike protein (S), 4 deletions and 6 synonymous. By the end of December, this variant had been detected in 31 countries and territories in 5 of the 6 WHO regions. One of the mutations is the deletion at positions 69-70 in the spike protein. Detection of the HV 69/70 deletion is of great importance since it has been related to immune leakage in immunosuppressed patients and to increased viral infectivity. Another cause for concern in relation to the HV 69/70 deletion is that it affects the sensitivity of virus detection using molecular techniques (RT-PCR) that detects the S gene.

The presence of the HV 69/70 deletion is associated with the UK variant, lineage B.1.1.7, however, other variants such as B.1.1.298 (Danish lineage) or B.1.258 also have this deletion.

The appearance of variants that increase the transmissibility of the virus, its virulence or that escape the action of the neutralizing antibodies generated after natural infection or the vaccine, constitute a first-order public health problem that can have an important impact on control of the pandemic. A concern regarding the new variants is that their detection by molecular techniques (RT-PCR) could be affected. For this reason, most commercial PCR tests use multiple targets located on different SARS-CoV-2 genes, so the likelihood that new variants will not be detected decreased.

### 3. Principle of the procedure

VIASURE SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Real Time PCR Detection Kit is designed for the qualitative detection of RNA from the SARS-CoV-2 and the HV 69/70 deletion of the S gene for SARS-CoV-2 in nasopharyngeal swabs. The detection is done in one step real time RT format where the reverse transcription and the subsequent amplification of specific target sequence occur in the same reaction well. The isolated RNA target is transcribed generating complementary DNA by reverse transcriptase which is followed by the amplification of a conserved region of the S gene for SARS-CoV-2 HV 69/70 deletion and by the amplification of a conserved region of the ORF1ab and N genes for SARS-CoV-2 using specific primers and a fluorescent-labeled probe.

VIASURE SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an

increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase and retrotranscriptase) in a stabilized format, as well as an **endogenous internal control** to monitor the extraction process and/or discard the inhibition of the polymerase activity. The assay uses a human housekeeping gene as an **endogenous Internal Control (IC)** (human RNase P gene). Human housekeeping genes are involved in basic cell maintenance and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels.

## 4. Reagents provided

VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit includes the materials and reagents detailed in Annex 1 for open format with internal control products and Annex 2 for tube format with internal control products.

## 5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- Real Time PCR compatible plastic consumables (i.e. individual tubes, well-strips and/or microplates). Only for Tubes format (Annex 2 and 4).
- RNA extraction kit.
- Collection and transport system.
- Laboratory freezers: - 30°C to - 10°C and/or ≤ -70°C.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.

VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-Time PCR System, Roche Molecular Diagnostics Cobas z480 Analyzer and NEOS-96 qPCR (Linear Chemicals). When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506). To check thermocycler compatibility and most common detection channels consult website [www.certest.es](http://www.certest.es).

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500\*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

\*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

## 6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from light.
- For Tube format kits: Once the SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Reaction-Mix tube has been reconstituted, it may be kept it at 25°C±5°C or 2°C to 8°C for up to 4 hours. For a longer period of time, it is recommended store at -20°C and to separate in aliquots to minimize freeze and thaw cycles (up to 6 times).

## 7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures (including training on the Real Time PCR instrument (thermocycler) and Nucleic acid extraction system).
- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (for references: VS-SUK213L, and VS-SUK213H). Remove any air excess in the pouches prior to closing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the functionality of the test.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.

- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE SARS-CoV-2 & UK Variant (*S* UK, ORF1ab and *N* genes) Real Time PCR Detection Kit and any additional reagents or equipment required for testing are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents at all times. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink, smoke or apply cosmetic products in the working area. Wash your hands after finishing the test.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- In accordance with Regulation (EC) No 1907/2006 (REACH), "VIASURE Real Time PCR Detection Kits" does not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment because it does not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

## 8. Test procedure

Please see Annex 1 for Open format with internal control products Test Procedure and Annex 2 for Tube format with internal control products Test Procedure.

### 8.1. Specimen collection, transport and storage

The VIASURE SARS-CoV-2 & UK Variant (*S* UK, ORF1ab and *N* genes) Real Time PCR Detection kit has been validated on nasopharyngeal specimens collected with synthetic fiber swabs with plastic and placed immediately into a sterile transport tube containing Universal transport medium (UTM). Other types of samples must be validated by the user.

Collection, storage and transport specimens should be maintained per the conditions validated by the user. Overall, respiratory samples should be collected and labelled appropriately in clean containers with or without transport media (depending on sample type) and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at 2 to 8°C for up to 72 hours, following the local and national regulations for the transport of pathogen material. For long term transport (more than 72 hours), It is recommended shipping at -20°C or lower. It is recommended to use fresh specimens for the test. The samples can

be stored at 2 to 8°C for up to 72 hours or frozen at -20°C or ideally at -70°C for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

The respiratory specimens must be collected, transport and storage according to appropriate laboratory guidelines. For details, refer to the CDC guideline (Specimen collection guidelines. Website <https://www.cdc.gov/urdo/downloads/SpecCollectionGuidelines.pdf>) and Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens for COVID-19. Website <https://www.cdc.gov/coronavirus/2019-ncov/lab/guidelines-clinical-specimens.html>) and the IDSA guideline (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94).

## 8.2. RNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For RNA extraction from respiratory samples, you can use your manual or automatic routine optimized system, or any commercially available RNA extraction kit and follow the manufacturer's instructions. The following extraction kits have been validated:

- MagMAX™ Viral/Pathogen II (MVP II) Nucleic Acid Isolation Kit using the KingFisher Flex System instrument (ThermoFisher).
- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).

## 9. Result interpretation

All the result of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Endogenous Internal Control (IC) signal to verify the extraction procedure and/or correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

**It is recommended to set the threshold values for each channel (target) independently by the end-user.** Use the Positive Control amplification curve as a starting point during the run validation (before than interpretation of patient sample results), in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal in the positive control well.

For a valid diagnostic test run, the following control conditions must be met:

Controls	S gene (HV 69/70 deletion) (FAM) <sup>1</sup>	ORF1ab gene (ROX) <sup>1</sup>	N gene (Cy5) <sup>1</sup>	Endogenous Internal Control (HEX)	Interpretation of Controls
<b>Positive Control (PC)</b>	≤40	≤40	≤40	≤40 <sup>2</sup>	<b>Valid</b>
<b>Negative Control (NC)</b>	≥40 or no signal	≥40 or no signal	≥40 or no signal	≥40 or no signal	<b>Valid</b>

Table 1. Expected Performance of Controls

1 In cases where one or more controls fail (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

2 The positive template control includes human housekeeping RNase P gene target; therefore, amplification signals are observed in all target channels, including the Endogenous Internal Control.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted.

For interpretation of patient sample results, use the following table, read and analyze the results:

<b>S gene (HV 69/70 deletion) (FAM)<sup>1</sup></b>	<b>ORF1ab gene (ROX)<sup>1</sup></b>	<b>N gene (Cy5)<sup>1</sup></b>	<b>Endogenous Internal Control (HEX)</b>	<b>Interpretation for patients' samples</b>	
≤40	≤40	≤40	≤40 or no signal <sup>1</sup>	<b>Valid</b>	<b>HV 69/70 deletion Detected</b>
≥40 or no signal	≤40	≤40	≤40 or no signal <sup>1</sup>	<b>Valid</b>	<b>SARS-CoV-2 Detected (HV 69/70 deletion not detected)</b>
≥40 or no signal	≥40 or no signal	≤40	≤40 or no signal <sup>1</sup>	<b>Inconclusive</b>	If only one SARS-CoV-2 target gene amplifies, repeat test depending on the available material:  a) obtain a new specimen, re-extract and retest (ideally) or,  b) re-extract and retest another aliquot of the same specimen or,  c) repeat RT-qPCR with the same isolated RNA sample.  After retesting one time, if at least one target gene is positive, the sample should be considered SARS-CoV-2 positive.
≥40 or no signal	≤40	≥40 or no signal	≤40 or no signal <sup>1</sup>		
≤40	≥40 or no signal	≥40 or no signal	≤40 or no signal <sup>1</sup>	<b>Inconclusive<sup>3</sup></b>	If S gene (HV 69/70 deletion) target gene amplifies, regardless of the amplification of other targets, repeat test depending on the available material:  a) obtain a new specimen, re-extract and retest (ideally) or,  b) re-extract and retest another aliquot of the same specimen or,  c) repeat RT-qPCR with the same isolated RNA sample.  After retesting one time, if S gene (HV 69/70 deletion) is positive, the sample should be considered HV 69/70 deletion detected.
≤40	≤40	≥40 or no signal	≤40 or no signal <sup>1</sup>		
≤40	≥40 or no signal	≤40	≤40 or no signal <sup>1</sup>		
≥40 or no signal	≥40 or no signal	≥40 or no signal	≤ 35 <sup>2</sup>	<b>Valid</b>	<b>Targets not Detected <sup>2</sup></b>
≥40 or no signal	≥40 or no signal	≥40 or no signal	≥ 35 or no signal <sup>2</sup>	<b>Invalid</b>	<b>Test Failure – Repeat Testing <sup>2</sup></b>

Table 2. Interpretation of individual patient sample results. Ct values. no signal = no amplification curve.

**1** The endogenous Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.

**2** In the case of SARS-CoV-2 target genes negative, IC must show an amplification signal with Ct less than 35. The Ct value could be very variable due to the Endogenous Internal Control is a human housekeeping gene that should be present in all human nucleated cells in the original sample. If there is an absence of signal or Ct value ≥ 35 of the endogenous Internal Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the RT-qPCR diluting the RNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.


**3** A sample only positive for the S gene target, can be considered as probable positive for SARS CoV-2, probably due to a low RNA load near the detection limit or to a different amplification efficiency between the targets. It is recommended to retest the sample to confirm the result (Consult Interpretation for patients' samples column in detail).

The presence of the HV 69/70 deletion is associated with the UK variant, lineage B.1.1.7, however, final assignment to a lineage must be done by sequencing.

In case of a continued ambiguous result, it is recommended to review the instructions for use; the extraction process used by the user; to verify the correct performance of each RT-qPCR steps and review the parameters; and to check the sigmoid shape of the curve and the intensity of fluorescence. It is also recommended to repeat the assay, preferably in duplicate. Depending on the available material:

- repeat RT-qPCR with the same isolated RNA sample, or
- re-extract and retest another aliquot of the same specimen or,
- obtain a new specimen and retest.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

## 10. Limitations of the test

- The results of the test should be evaluated by a healthcare professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with RNA extracted from respiratory samples (nasopharyngeal swab).
- The quality of the test depends on the quality of the sample; properly extracted nucleic acid from clinical samples must be extracted.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by SARS-CoV-2, either the high number of cDNA template copies which contains each SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Positive Control vial, samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.
- The specific primer and probe combinations for detection of the ORF1ab and N genes used in VIASURE SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Real Time PCR Detection Kit designed for the

detection of SARS-CoV-2, do not show significant combined homologies with the human genome, human microflora, or other coronaviruses (with the exception of some N and/or ORF1ab sequences from SARS-CoV, and other coronaviruses identified in bats and pangolin), which might result in predictable false positive.

- False Negative results may arise from several factors and their combinations, including:
  - Improper specimens' collection, transport, storage, and/or handling methods.
  - Improper processing procedures (including RNA extraction).
  - Degradation of the viral RNA during sample shipping/storage and/or processing.
  - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown SARS-CoV-2 variants.
  - A viral load in the specimen below the limit of detection for the assay.
  - The presence of RT-qPCR inhibitors or other types of interfering substances. The impacts of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent COVID-19 or used during the treatment of the infection have not been evaluated.
  - Failure to follow instructions for use and the assay procedure.

If in doubt, refer to section 9 to check the correct interpretation of the results.

- A positive test result does not necessarily indicate the presence of viable viruses and does not imply that these viruses are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of targets viral sequences.
- The presence of the HV 69/70 deletion is associated with the UK variant, lineage B.1.1.7, however, final assignment to a lineage must be done by sequencing.
- Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak viral levels during infections caused by SARS-CoV-2 have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the virus.
- If diagnostic tests for other respiratory illnesses are negative and the patient's clinical presentation and epidemiological information suggest that SARS-CoV-2 infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.
- Some samples may fail to exhibit RNase P amplification curves due to low human cell numbers in the original clinical sample. A negative IC signal does not preclude the presence SARS-CoV-2 RNA in a clinical specimen.

## 11. Quality control

VIASURE SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the endogenous internal control (IC) in each well confirms the correct performance of the technique.

## 12. Performance characteristics

### 12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit was tested using respiratory clinical samples (nasopharyngeal swabs) from patients with suspected respiratory infection. The results were as follows:

	<b>Site</b>	<b>Sample type</b>	<b>Workflow</b>	<b>Target</b>
1	Fundación Jiménez Díaz (Madrid, Spain)	nasopharyngeal swab	MagMAX™ Viral/Pathogen II (MVP II) Nucleic Acid Isolation Kit using the KingFisher Flex System instrument (ThermoFisher) + Bio-Rad CFX96™ Real-Time PCR Detection System	HV 69/70 deletion
				SARS-CoV-2 (non HV 69/70 deletion)
2	"Servicio de Microbiología" of the Hospital Universitario Marqués de Valdecilla (Santander, Spain)	nasopharyngeal swab	MagMAX™ Viral/Pathogen II (MVP II) Nucleic Acid Isolation Kit using the KingFisher Flex System instrument (ThermoFisher) + Bio-Rad CFX96™ Real-Time PCR Detection System	HV 69/70 deletion
				SARS-CoV-2 (non HV 69/70 deletion)

Table 3. Site, sample type, workflow and target.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE SARS-CoV-2 & UK Variant (*S* UK, *S* and *N* genes) Real Time PCR Detection Kit were calculated in relation to each comparator assay as shown in the following tables:

<b>Site</b>	<b>Comparator assay</b>	<b>Target</b>	<b>TP</b>	<b>TN</b>	<b>FP</b>	<b>FN</b>	<b>Sensitivity</b>	<b>Specificity</b>	<b>PPV</b>	<b>NPV</b>
1	TaqPath COVID-19 CE-IVD RT-PCR Kit molecular assay + sequencing	HV 69/70 deletion	22	38	0	0	1 (0.84-1)	1 (0.90-1)	1 (0.84-1)	1 (0.90-1)
		SARS-CoV-2 (non HV 69/70 deletion)	18	42	0	0	1 (0.78-1)	1 (0.91-1)	1 (0.78-1)	1 (0.91-1)
2	TaqPath COVID-19 CE-IVD RT-PCR Kit molecular assay + sequencing	HV 69/70 deletion	50	49	0	0	1 (0.93-1)	1 (0.92-1)	1 (0.93-1)	1 (0.92-1)
		SARS-CoV-2 (non HV 69/70 deletion)	49	50	0	0	1 (0.92-1)	1 (0.93-1)	1 (0.92-1)	1 (0.93-1)

Table 4. True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV values for VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit.

Results show high agreement to detect SARS-CoV-2 and the HV 69/70 deletion using VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit.

## 12.2. Analytical sensitivity

VIASURE SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Real Time PCR Detection Kit has a detection limit (LoD) of 40 genome copies/rxn for S gene (HV 69/70 deletion), 40 genome copies/rxn for ORF1ab gene and 80 genome copies/rxn for N gene. (Figures 2, 3 and 4).

Figure 2. Dilution series of S gene (HV 69/70 deletion) ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (FAM channel).


Figure 3. Dilution series of ORF1ab gene ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (ROX channel).



Figure 4. Dilution series of N gene ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Cy5 channel).

## 12.3. Analytical specificity

The specificity of the SARS-CoV-2 assay was confirmed by testing a panel consisting of different microorganisms representing the most common respiratory pathogens. No cross-reactivity was detected between any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing					
Human Adenovirus types 1-5, 8, 15, 31, 40 and 41	-	Influenza A/New Caledonia/20/99(H1N1) virus	-	<i>Legionella longbeachae</i>	-
Human Bocavirus	-	Influenza A/California/7/2009(H1N1)pdm09	-	<i>Legionella micdadei</i>	-
<i>Bordetella bronchiseptica</i>	-	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	<i>Legionella pneumophila</i>	-
<i>Bordetella holmesii</i>	-	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	-	Human metapneumovirus A and B	-
<i>Bordetella parapertussis</i>	-	Influenza A/Victoria/210/2009 (H3N2)	-	<i>Moraxella catarrhalis</i>	-
<i>Bordetella pertussis</i>	-	Influenza A/Thüringen/5/17 (H3N2) virus	-	<i>Mycoplasma pneumoniae</i>	-
<i>Chlamydia caviae</i>	-	Influenza A/Switzerland/9715293/2013 (H3N2) virus	-	<i>Mycobacterium tuberculosis</i> not rifampin resistant	-
<i>Chlamydia psittaci</i> genotype A and C	-	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	-	Human parainfluenza 1, 2, 3 and 4 viruses	-
<i>Chlamydophila pneumoniae</i> CM-1	-	Influenza A/South Australia/55/2014, IVR-175 (H3N2) virus	-	<i>Pneumocytis jirovecii</i>	-
Human coronavirus 229E, OC43, NL63 and HKU1	-	Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus	-	Human rhinovirus type C	-
MERS Coronavirus	-	Influenza A/Anhui/1/2013 (H7N9) virus	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
SARS Coronavirus Strain Frankfurt 1	-	Influenza B/Brisbane/60/2008	-	<i>Staphylococcus epidermidis</i>	-
Enterovirus 68 and 71	-	Influenza B/Florida/04/06 virus	-	<i>Streptococcus pneumoniae</i> Z202	-
Enterovirus Echovirus 11 and 30	-	Influenza B/Phuket/3073/2013 virus	-	<i>Streptococcus pyogenes</i>	-
Enterovirus Coxsackievirus A24, A9 and B3	-	<i>Legionella bozemanii</i>	-	<i>Streptococcus salivarius</i>	-
<i>Haemophilus influenzae</i> MinnA	-	<i>Legionella dumoffii</i>	-	Respiratory syncytial virus (RSV) A and B	-

Table 5. Reference pathogenic microorganisms used in this study

## 12.4. Analytical reactivity

The reactivity of the VIASURE SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Real Time PCR Detection Kit for was evaluated against RNA from Human 2019-nCoV strain BetaCoV/Germany/BavPat1/2020 p.1, Human 2019-nCoV strain 2019-nCoV/Italy-INMI1, SARS-CoV-2 strain 2019nCoV/USA-WA1/2020, SARS-CoV-2 BetaCoV/Berlin/ChVir1670/2020\_IsolatBER, SARS-CoV-2 BetaCoV/Munich/ChVir984/2020, SARS-CoV-2 BetaCoV/Baden-Wuerttemberg/1/ChVir1577/2020\_IsolatBER and synthetic RNA controls for four variants of the SARS-CoV-2 virus: MT007544.1 (SARS-CoV2 isolate Australia/VIC01/2020), MN908947.3 (SARS-CoV-2 isolate Wuhan-Hu-1), B.1.1.7\_710528 (Twist Synthetic SARS-CoV-2 RNA Control 14) and B.1.1.7\_601443 (Twist Synthetic SARS-CoV-2 RNA Control 15), showing positive result.

## ANNEX 1

### OPEN FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE SARS-CoV-2 & UK Variant ( <i>S</i> UK, <i>ORF1ab</i> and <i>N</i> genes) Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-SUK206L
VIASURE SARS-CoV-2 & UK Variant ( <i>S</i> UK, <i>ORF1ab</i> and <i>N</i> genes) Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-SUK206H
VIASURE SARS-CoV-2 & UK Variant ( <i>S</i> UK, <i>ORF1ab</i> and <i>N</i> genes) Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-SUK212L
VIASURE SARS-CoV-2 & UK Variant ( <i>S</i> UK, <i>ORF1ab</i> and <i>N</i> genes) Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-SUK212H
VIASURE SARS-CoV-2 & UK Variant ( <i>S</i> UK, <i>ORF1ab</i> and <i>N</i> genes) Real Time PCR Detection Kit 96-well plate, low profile	VS-SUK213L
VIASURE SARS-CoV-2 & UK Variant ( <i>S</i> UK, <i>ORF1ab</i> and <i>N</i> genes) Real Time PCR Detection Kit 96-well plate, high profile	VS-SUK213H
VIASURE SARS-CoV-2 & UK Variant ( <i>S</i> UK, <i>ORF1ab</i> and <i>N</i> genes) Real Time PCR Detection Kit 1 x 8-well strips, low profile	VS-SUK201L
VIASURE SARS-CoV-2 & UK Variant ( <i>S</i> UK, <i>ORF1ab</i> and <i>N</i> genes) Real Time PCR Detection Kit 1 x 8-well strips, high profile	VS-SUK201H

Table A1 1. References

### A1.1 Principle of the procedure

VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase and retrotranscriptase) in a stabilized format, as well as an endogenous internal control to monitor the extraction process and/or discard the inhibition of the polymerase activity. The assay uses a human housekeeping gene as an endogenous Internal Control (IC) (human RNase P gene). Human housekeeping genes are involved in basic cell maintenance and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels.

Target	Channel	Gene
SARS-CoV-2 (HV 69/70 deletion)	FAM	<i>S</i> gene (HV 69/70 deletion)
SARS-CoV-2 ( <i>ORF1ab</i> gene)	ROX	<i>ORF1ab</i> gene
SARS-CoV-2 ( <i>N</i> gene)	Cy5	<i>N</i> gene
Endogenous Internal control (IC)	HEX, VIC or JOE *	human RNase P gene

Table A1 2. Target, channel and genes.

\*Depending on the equipment used select the proper detection channel, to check most common detection channels consult the website [www.certest.es](http://www.certest.es).

### A1.2 Reagents provided

VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables A1.3 and A1.4. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table A1.3 includes materials and reagents to be used with 8-well strips compatible

devices. Table A1.4 includes materials and reagents to be used with 96-well plate compatible devices. (Consult the thermocycler compatibility on CerTest's website [www.certest.es](http://www.certest.es)).

Reagent/Material	Description	Colour	Amount
SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous Internal control in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNAse free	RNAse/DNAse free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 x 8-cap strip

Table A1 3. Reagents and materials provided in VIASURE SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Real Time PCR Detection Kit with Ref. VS-SUK206L, VS-SUK206H, VS-SUK212L and VS-SUK212H.

Reagent/Material	Description	Color	Amount
SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNAse free	RNAse/DNAse free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table A1 4. Reagents and materials provided in VIASURE SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Real Time PCR Detection Kit with Ref VS-SUK213L and VS-SUK213H.

## A1.3 Test procedure

### A1.3.1 Lyophilized positive control

SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNAse free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

### A1.3.2 PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of RNA sample, reconstituted SARS-CoV-2 & UK Variant (S UK, ORF1ab and N genes) Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate.

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website [www.certest.es](http://www.certest.es)).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Reverse transcription	15 min	45°C
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table A1 5. PCR protocol

Fluorogenic data should be collected during the extension step (\*) through the FAM (HV 69/70 deletion), ROX (ORF1ab gene), Cy5 (N gene) and HEX, JOE or VIC (Endogenous Internal Control (IC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website [www.certest.es](http://www.certest.es)). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005PTM Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

## ANNEX 2

### TUBE FORMAT WITH INTERNAL CONTROL

Annex for the following references:

PRODUCT	REFERENCE
VIASURE SARS-CoV-2 & UK Variant ( <i>S</i> UK, <i>ORF1ab</i> and <i>N</i> genes) Real Time PCR Detection Kit, 4 tubes x 24 reactions	VS-SUK296T

Table A2. 1. References.

### A2.1 Principle of the procedure

VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit contains in each Reaction-Mix tube all the components necessary for 24 real time PCR reactions (specific primers/probes, dNTPS, buffer, polymerase and retrotranscriptase) in a stabilized format, as well as an endogenous internal control to monitor the extraction process and/or discard the inhibition of the polymerase activity. The assay uses a human housekeeping gene as an endogenous Internal Control (IC) (human RNase P gene). Human housekeeping genes are involved in basic cell maintenance and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels.

Target	Channel	Gene
SARS-CoV-2 (HV 69/70 deletion)	FAM	<i>S</i> gene (HV 69/70 deletion)
SARS-CoV-2 ( <i>ORF1ab</i> gene)	ROX	<i>ORF1ab</i> gene
SARS-CoV-2 ( <i>N</i> gene)	Cy5	<i>N</i> gene
Endogenous Internal control (IC)	HEX, VIC or JOE *	human RNase P gene

Table A2. 2.Target, channel and genes.

\*Depending on the equipment used select the proper detection channel, channel, to check most common detection channels consult website [www.certest.es](http://www.certest.es).

### A2.2 Reagents provided

VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit includes the following materials and reagents detailed in Table A2.3.

Reagent/Material	Description	Colour	Amount
SARS-CoV-2 & UK Variant ( <i>S</i> UK, <i>ORF1ab</i> and <i>N</i> genes) Reaction-Mix tube	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous Internal control in stabilized format	White	4 vials
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
SARS-CoV-2 & UK Variant ( <i>S</i> UK, <i>ORF1ab</i> and <i>N</i> genes) Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNAse free	RNAse/DNAse free water	White	1 vial x 1 mL

Table A2. 3. Reagents and materials provided in VIASURE SARS-CoV-2 & UK Variant (*S* UK, *ORF1ab* and *N* genes) Real Time PCR Detection Kit with Ref. VS-SUK296T.

## A2.3 Test procedure

### A2.3.1 Lyophilized positive control

SARS-CoV-2 & UK Variant (*S* UK, ORF1ab and *N* genes) Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized SARS-CoV-2 & UK Variant (*S* UK, ORF1ab and *N* genes) Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNAse free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. It is recommended to separate it in aliquots to minimize freeze and thaw cycles.

### A2.3.2 Lyophilized reaction mix tube

Determine the number of required reactions including samples and controls (one positive and negative control must be included in each run). Obtain the correct number of lyophilized Reaction-Mix vials (24-reactions each one) for testing.

Recommendation is to open and manipulate the SARS-CoV-2 & UK Variant (*S* UK, ORF1ab and *N* genes) Reaction-Mix tube in pre-PCR laboratory area. Open lyophilized Reaction-mix tube (white vial) carefully to avoid disruption of the pellet and add 390 µL of Rehydration Buffer (blue vial) supplied. Mix gently by pipetting up and down. Spin down briefly to remove bubbles generated during mixing.

Once the Reaction-Mix tube has been re-suspended, return unused reagents to the appropriate storage conditions at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The volume of the rehydrated Reaction-Mix is sufficient for 24 reactions. The rehydrated Reaction-Mix may be kept at 25°C±5°C or 2-8°C for up to 4-hours (see Transport and storage conditions section for additional storage options).

### A2.3.3 PCR protocol

- 1) Adding rehydrated Reaction-Mix to the number of required wells.

Add 15 µL of rehydrated SARS-CoV-2 & UK Variant (*S* UK, ORF1ab and *N* genes) Reaction-Mix (white vial) into each tube.

- 2) Adding samples and controls.

Add 5 µL of RNA sample, reconstituted SARS-CoV-2 & UK Variant (*S* UK, ORF1ab and *N* genes) Positive Control (red vial) or Negative Control (violet vial) in different wells and close the tubes with caps or seal the plate. Centrifuge briefly.

Load the plate, the strips, or tube in the thermocycler.

- 3) Set up the thermocycler (consult thermocycler compatibility on CerTest's website [www.certest.es](http://www.certest.es)).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Reverse transcription	15 min	45°C
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table A2. 4. PCR protocol.

Fluorogenic data should be collected during the extension step (\*) through the FAM (HV 69/70 deletion), ROX (ORF1ab gene), Cy5 (N gene) and HEX, JOE or VIC (Endogenous Internal Control (IC)). Depending on the equipment used select the proper detection channel (to check most common detection channels consult website [www.certest.es](http://www.certest.es)). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

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## Bibliography/Bibliografía

1. Huang, C. et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. *The Lancet*, 2020. DOI : 10.1016/S0140-6736(20)30183-5.
2. Zhu N. et al. A novel coronavirus from patients with pneumonia in China, 2019. *New England Journal of Medicine.*, 2020. DOI : 10.1056/NEJMoa2001017.
3. World Health Organization. MERS situation update. January 2020. Available from <https://applications.emro.who.int/docs/EMCSR254E.pdf?ua=1> Accessed January 2021.
4. Chen N. et al.. Epidemiological and Clinical Characteristics of 99 Cases of 2019-Novel Coronavirus (2019-nCoV) Pneumonia in Wuhan, China. *The Lancet*, 2020. DOI: 10.1016/S0140-6736(20)30211-7.
5. Lv D.F. et al. Dynamic change process of target genes by RT-PCR testing of SARS-CoV-2 during the course of a Coronavirus Disease 2019 patient. *Clinica Chimica Acta* 2020; 506: 172-175.
6. World Health Organization. Clinical management of COVID-19 disease" Interim guidance 27 May 2020. Available from <https://www.who.int/publications/item/clinical-management-of-covid-19> Accessed January 2021.
7. Lu R. et al. Genomic characterisation and epidemiology of 2019 novel coronavirus: implications for virus origins and receptor binding. *The Lancet*, 2020. DOI : 10.1016/S0140-6736(20)30251-8.
8. Rothe C. et al. Transmission of 2019-nCoV Infection from an Asymptomatic Contact in Germany. *New England Journal of Medicine*, 2020. DOI : 10.1056/NEJMc2001468.
9. Centers of Disease Control and Prevention (CDC). Coronavirus Disease 2019 (COVID-19), Symptoms of Coronavirus. Available from <https://www.cdc.gov/coronavirus/2019-ncov/about/symptoms.html> Accessed January 2021.
10. Centers of Disease Control and Prevention (CDC). Coronavirus Disease 2019 (COVID-19), Older Adults. Available from <https://www.cdc.gov/coronavirus/2019-ncov/need-extra-precautions/older-adults.html> Accessed January 2021.
11. World Health Organization. Laboratory testing for 2019 novel coronavirus (2019-nCoV) in suspected human cases. Interim guidance. 19 March 2020. Available from <https://www.who.int/publications-detail/laboratory-testing-for-2019-novel-coronavirus-in-suspected-human-cases-20200117> Accessed January 2021.
12. Yan Y et al. Laboratory testing of SARS-CoV, MERS-CoV, and SARS-CoV-2 (2019-nCoV): Current status, challenges, and countermeasures. *Reviews in Medical Virology* 2020; 30(3):e2106.
13. Centers of Disease Control and Prevention (CDC). 2019-Novel Coronavirus (2019-nCoV) Real-time rRT-PCR Panel Primers and Probes. <https://www.cdc.gov/coronavirus/2019-ncov/downloads/rt-pcr-panel-primer-probes.pdf> Accessed January 2021.
14. Chu D.K.W. et al. Molecular Diagnosis of a Novel Coronavirus (2019-nCoV) Causing an Outbreak of Pneumonia. *Clinical Chemistry* 2020;66(4): 549-555.
15. Corman V.M. et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. *European communicable disease bulletin* 2020;25(3).
16. European Centre for Disease Prevention and Control. Novel coronavirus disease 2019 (COVID-19) pandemic: increased transmission in the EU/EEA and the UK – sixth update – 12 March 2020. Stockholm: ECDC; 2020. Available from <https://www.ecdc.europa.eu/sites/default/files/documents/RRA-sixth-update-Outbreak-of-novel-coronavirus-disease-2019-COVID-19.pdf> Accessed January 2021.
17. Lim, Y. X., Ng, Y. L., Tam, J. P., & Liu, D. X. (2016). Human coronaviruses: a review of virus–host interactions. *Diseases*, 4(3), 26.

18. McBride R. et al. The coronavirus nucleocapsid is a multifunctional protein. *Viruses* 2014; 6(8):2991-3018.
19. Sheikh A. et al. Analysis of preferred codon usage in the coronavirus N genes and their implications for genome evolution and vaccine design. *Journal of Virological Methods* 2020; 277:113806.
20. World Health Organization. Public health surveillance for COVID-19. 16 December 2020. Available from [https://www.who.int/publications-detail/global-surveillance-for-human-infection-with-novel-coronavirus-\(2019-ncov\)](https://www.who.int/publications-detail/global-surveillance-for-human-infection-with-novel-coronavirus-(2019-ncov)) Accessed January 2021.
21. G. Neumann et al. Transmission of Influenza A viruses. *Virology* 2015; 234-246.
22. Y. Yang et al. Simultaneous typing and HA/NA subtyping of influenza A and B viruses including the pandemic influenza A/H1N1 2009 by multiplex real-time RT-PCR. *Journal of Virological Methods* 2010; 167(1): 37-44.
23. R.L. Kuo et al. Influenza A/B virus detection and influenza A virus subtyping with emphasis on the novel H7N9 virus by using multiplex real-time RT-PCR. *Journal of Virological Methods* 2014; 208:41-46.
24. W.P. Glezen et al. The burden of Influenza B: A structured literature review. *American Journal of Public Health* 2013; 103(3): 43-51.
25. F. de-Paris et al. Optimization of one-step duplex real-time RT-PCR for detection of influenza and respiratory syncytial virus in nasopharyngeal aspirates. *Journal of Virological Methods* 2012; 186(1-2): 189- 192.
26. M. Hindiyeh et al. Evaluation of Simplexa Flu A/B & RSV for direct detection of influenza viruses (A and B) and respiratory syncytial virus in patient respiratory samples. *Journal of Clinical Microbiology* 2013; 51(7): 2421-2424 .
27. World Health Organization. WHO information for molecular diagnosis of influenza virus—update. Available: [http://www.who.int/influenza/gisrs\\_laboratory/molecular\\_diagnosis/en/](http://www.who.int/influenza/gisrs_laboratory/molecular_diagnosis/en/). Accessed 2015 Dec 30.
28. Enfermedad por coronavirus, COVID-19, Información Científica-técnica. Centro de Coordinación de Alertas y Emergencias Sanitarias. Ministerio de Sanidad, España. 01-2021.

## Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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## Trademarks

- CFX™ and IQ5™ are trademarks of Bio-Rad Laboratories.
- COBAS® are registered trademark of Roche.
- AriaMx are trademarks of Agilent Technologies.

- MagMAX™ and KingFisher™ are trademarks property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.
- MagDEA® Dx and magLEAD® are registered trademark of PSS.

All other trademarks that may appear in this package insert are the property of their respective owners.

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**bestbion dx GmbH**