
COVID-19 IgM/IgG Combo FIA

STANDARD™ F COVID-19 IgM/IgG Combo FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

KIT CONTENTS

SPECIMEN COLLECTION AND PREPARATION

- Serum**
1. Collect the whole blood into the commercially available plain tube, NOT containing anti-coagulants such as heparin, EDTA or sodium citrate, by venipuncture and leave to settle for 30 minutes for blood clotting to occur.
2. If serum in the plain tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 1 week after collection. Using the specimen in the long-term keeping more than 1 week can cause non-specific reaction. For prolonged storage, it should be stored at below -40°C/-40°F. They should be brought to room temperature prior to use.

- Plasma**
1. Collect the venous blood into the commercially available EDTA anti-coagulant tube by venipuncture and centrifuge blood to get plasma specimen.
2. If plasma in the EDTA tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 1 week after collection. Using the specimen in the long-term keeping more than 1 week can cause non-specific reaction. For prolonged storage, it should be stored at below -40°C/-40°F. They should be brought to room temperature prior to use.

Whole blood

- Capillary whole blood**
1. Capillary whole blood should be collected aseptically by fingertip.
2. Clean the area to be pierced with an alcohol swab.
3. Squeeze the end of the fingertip and pierce with a sterile lancet.
4. Using a capillary tube, collect the 20μl of capillary whole blood to the black line of the capillary tube.
Venous whole blood
1. Collect the venous whole blood into the commercially available anti-coagulant tube such as EDTA by following the instructions.
2. If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 1 day after collection.
3. Do not use hemolyzed blood specimens.


Use separate disposable materials for each specimen in order to avoid cross-contamination which can cause erroneous results.

PREPARATION AND TEST PROCEDURE

- Preparation**
1. Load test device and collected specimen to room temperature (15-30°C / 59-86°F) prior to testing.
2. Carefully read instructions for use.
3. Read the instructions for use.

* Check the valid expiry date at the back of the foil pouch. Do not use if the expiry date has passed.

* The valid color of indicator is yellow.

The violet-line on the membrane of unused test device will disappear after use.

Do not write on the barcode or damage the barcode of the test device.

Analysis of specimen

- 'STANDARD' mode**
1. Prepare a STANDARD F Analyzer and select the "Standard Mode" according to the analyzer's manual.

STANDARD F200 analyzer *Workplace* → Run Test → Insert patient ID and / or operator ID → Standard Test mode → Insert patient ID and / or operator ID → On the analyzer

- 2. Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of the test device.

- 3. Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer will read the barcode data, and check the test device is valid.

- 4. Dispense the collected serum and plasma (10μl) or whole blood (20μl) to the specimen well of the test device.

- 5. Apply 3 drops of buffer to the specimen well of the test device.

- 6. After applying the buffer, immediately press the "TEST START" button.

- 7. The analyzer will automatically display the test result in 15 minutes.

'READ ONLY' mode

STANDARD™ F100 and F200 analyzer

- Take the test device out of the foil pouch and place it on a flat and dry surface. Write information on the back of the test device.
- Dispense the collected serum and plasma (10μl) or whole blood (20μl) to the specimen well of the test device.
- Apply 3 drops of buffer to the specimen well of the test device.
- Incubate the test device for 15 minutes outside of the analyzer.
- Prepare a STANDARD F Analyzer and select the "Read Only" mode according to the analyzer's manual.
- Insert the test device to the test slot of the analyzer. The analyzer automatically reads the barcode information of the test device and outputs the resulting values based on the test information.

QUALITY CONTROL

■ STANDARD F Analyzers calibration check

The calibration set test of STANDARD F Analyzers should be conducted according to the analyzer's manual.

When to use calibration set

- Before using for the first time
- When changing the analyzer
- After cleaning the test slot with your result
- If you want to check the performance of an analyzer and test device

How to use calibration set

Calibration set test is a mandatory function that ensures optimal performance by checking the internal components of the test device.

1. Select the "Calibration" menu.

2. The specific calibration set is included with the analyzer.
3. Insert the CAL-1 first, and then insert the CAL-2 for UV-LED testing and the CAL-3 for RGB-LED testing in order.

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BIBLIOGRAPHY

1. Clinical management of severe acute respiratory infection when novel coronavirus(nCoV) infection is suspected. Interim guidance. WHO.2020
2. Diagnostic detection of Wuhan coronavirus 2019 by real-time RT-PCR.2020
3. Diagnosis and treatment of pneumonia caused by new coronaviruses (trial version 4) National Health Commission. 2020
4. Guo L et al. Profiling Early Humoral Response to Novel Coronavirus Disease (COVID-19). Clinical Infectious Disease. 2020
5. Zhao J et al. Antibody responses to SARS-CoV-2 in patients of novel coronavirus disease 2019. Clinical Infectious Disease. 2020

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