
EN Q-NCOV-01G Cat.No.: 09COV30D


STANDARD Q COVID-19 Ag

STANDARD™ Q COVID-19 Ag Test

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

KIT CONTENTS**PREPARATION AND TEST PROCEDURE****■ PREPARATION**

1. Carefully read instructions for using the STANDARD Q COVID-19 Ag Test.
2. Check the expiry date at the back of the foil pouch. Do not use the kit, if expiry date has passed.
3. Check the test device and the desiccant pack in the foil pouch.

Foil pouch

Test device

Desiccant

■ COLLECTION OF SPECIMEN [Nasopharyngeal swab]

1. Insert a sterile swab into the nostril of the patient so that the swab reaches the depth equal to distance from nostrils to outer opening of the ear. Rotate the swab 3-4 times against the surface of the nasopharyngeal and withdraw the swab from the nasal cavity.

2. Insert the swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab more than 5 times.

3. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.

4. Press the nozzle cap tightly onto the tube.

5. Place the test device on a flat surface.

6. Dispense the specimen at 90 degree angle to allow for free falling drops and avoid bubbles.

7. Do not read test results after 30 minutes. It may give false results.

■ INTERPRETATION OF TEST RESULT*** "C" Control Line | I "T" Test Line**

Negative

Positive

Invalid

1. A purple colored band will appear in the top section of the result window to show that the test is working properly. This band is control line (C).

2. A purple colored band will appear in the lower section of the result window. This band is test line of SARS-CoV-2 antigen (T).

3. Even if the control line is faint, or the test line isn't uniform, the test should be considered to be performed properly and the test result should be interpreted as a positive result.

* The presence of any line no matter how faint the result is considered positive.

* Positive results should be considered in conjunction with the clinical history and other data available.

■ EXPLANATION AND SUMMARY**■ Introduction**

Coronavirus is a single-stranded positive-sense virus with an envelope of about 80 to 120 nm in diameter. Its genetic material is the largest of all RNA viruses and is an important pathogen of many domestic animals, pets, and human diseases. It can cause a variety of acute and chronic diseases. Common signs of a person infected with a coronavirus include respiratory symptoms, fever, cough, shortness of breath, and dyspnea. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure, and even

death.

The new 2019 coronavirus, or "SARS-CoV-2 (COVID-19)", was discovered because of Wuhan Viral Pneumonia cases in 2019, and was named by the World Health Organization on January 12, 2020, confirming that it can cause colds and the Middle East Respiratory Syndrome (MERS) and more serious diseases such as acute respiratory syndrome (SARS). This kit is helpful for the auxiliary diagnosis of coronavirus infection. The test results are for clinical reference only and cannot be used as a basis for confirming or excluding cases alone.

■ Standard use

STANDARD Q COVID-19 Ag Test is a rapid chromatographic immunoassay for the qualitative detection of specific antigens to SARS-CoV-2 present in human nasopharyngeal specimens. This product is intended for healthcare professionals at the clinical setup and point of care sites, as an aid to early diagnosis of SARS-CoV-2 infection in patient with clinical symptoms of SARS-CoV-2 infection. It provides only an initial screening test result. This product is strictly for medical professional use only and not intended for personal use. The administration of the test and the interpretation of the results should be done by a trained health professional. The result of this test should not be the sole basis for the diagnosis; confirmatory testing is required.

■ Principle

STANDARD Q COVID-19 Ag Test has two pre-coated lines, "C" Control line, "T" Test line on the surface of the nitrocellulose membrane. Both the control line and test line in the result window are not visible before applying any specimens. Mouse monoclonal anti-SARS-CoV-2 antibody is coated on the test line region and mouse monoclonal anti-Chicken IgY antibody is coated on the control line region. Mouse monoclonal anti-SARS-CoV-2 antibody conjugated with color particles are used to detect SARS-CoV-2 antigen in the specimen. If SARS-CoV-2 antigen is present in the specimen, then color appears in the test line. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

7. Positive test results do not rule out co-infections with other pathogens.

8. When using VTM, sensitivity can be reduced due to dilution.

9. Only Copan UTM, BD UTU and STANDARD™ Transport Medium have been validated with the assay.

■ KIT STORAGE AND STABILITY

Store the kit at 2-30°C / 36-86°F out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

■ WARNINGS AND PRECAUTIONS

1. Bring the kit contents and the specimen to room temperature before testing.

2. Do not re-use the test kit.

3. Do not use the test kit if the pouch is damaged or the seal is broken.

4. Do not use the extraction buffer tube of another lot.

5. Do not smoke, drink or eat while handling specimen.

6. Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.

7. Clean up spills thoroughly using an appropriate disinfectant.

8. Handle all specimens as if they contain infectious agents.

9. Observe established precautions against microbiological hazards throughout

10. Dispose of all specimens and materials used to perform the test as biohazard waste. Laboratory chemical and biohazard wastes must be handled and discarded in accordance with all local, state, and national regulations.

11. Desiccant in foil pouch is to absorb moisture and keep humidity from

the test kit contents.

12. When using viral transport medium (VTM), it is important to ensure that the VTM containing the specimen is warmed to room temperature. Cold specimens will not flow correctly and can lead to erroneous or invalid results. Several minutes will be required to bring a cold specimen to room temperature.

PERFORMANCE CHARACTERISTICS**[Clinical evaluation]**

The prospective diagnostic evaluation of STANDARD Q COVID-19 Ag Test with a total number of enrolled individuals of 1659 was conducted by FIND with collaborators in Germany and Brazil.

A total of 153 positive specimens from Germany and Brazil were tested using the STANDARD Q COVID-19 Ag Test. These specimens consisted of nasopharyngeal swabs from symptomatic patients. The specificity of STANDARD Q COVID-19 Ag Test was tested using 1506 negative specimens.

The sensitivity and specificity of the STANDARD Q COVID-19 Ag Test was compared to the site-specific RT-PCR method. The pooled sensitivity was 84.97% (130/153, 95% CI 78.3% - 90.23%) and the pooled specificity was 98.94% (1490/1506, 95% CI 98.28% - 99.39%). Performance data was calculated from a study of patients within 24 days of onset of symptoms.

Table 1. STANDARD Q COVID-19 Ag Test result by FIND.

Country	Brazil	Germany	Overall
Sensitivity (Ct < 25)	95.92% (47/49, 95% CI 86.02-99.50%)	100% (21/21, 95% CI 83.89-100%)	97.14% (68/70, 95% CI 90.06-99.65%)
Sensitivity (Ct < 33)	91.92% (91/99, 95% CI 84.70-96.45%)	87.80% (36/41, 95% CI 73.80-95.92%)	90.71% (127/140, 95% CI 84.64-94.96%)
(0 ≤ from the symptom onset days ≤ 3)	95% (19/20, 95% CI 75.13-99.87%)	85.71% (18/21, 95% CI 63.66-96.55%)	90.24% (37/41, 95% CI 76.87 - 97.28%)
Sensitivity (from the symptom onset days ≤ 7)	90.7% (88/97, 95% CI 83.12-95.67%)	80% (28/35, 95% CI 63.06-91.56%)	87.88% (116/132, 95% CI 81.06-92.91%)
Clinical Sensitivity	88.68%	76.6%	84.97%
Clinical Specificity	97.6% (94/106, 95% CI 81.06-94.01%)	99.3% (36/47, 95% CI 61.97-87.70%)	98.94% (130/153, 95% CI 78.3 - 90.23%)
Clinical Specificity	97.6% (287/294, 95% CI 95.2-98.8%)	99.3% (1203/1212, 95% CI 98.6-99.6%)	98.94% (1490/1506, 95% CI 98.28-99.39%)

ANALYTICAL PERFORMANCE**1. Limit of Detection (LoD)**The SARS-CoV-2 positive specimen was prepared by spiking Inactivated SARS-CoV-2 (2019-nCoV) NCCP 43326/2020/Korea strain to SARS-CoV-2 negative nasopharyngeal swab confirmed with PCR. LoD is determined as $3.12 \times 10^{2.2} \text{ TCID}_{50}/\text{ml}$ for direct Nasopharyngeal swab, $5 \times 10^{2.2} \text{ TCID}_{50}/\text{ml}$ for Nasopharyngeal swab stored in VTM by testing serially diluted the mock positive specimen.**2. Cross-Reactivity+ Microbial interference**

There was no cross-reaction and interference with the potential cross-reacting microorganisms listed below except SARS-CoV.

Potential cross reacting substance	Strain	Concentration of potentially cross reacting substance
SARS-coronavirus	Urbani	3.5 pg/ml
MERS-Coronavirus	Florida/USA-2_Saudi Arabia_2014	$4 \times 10^4 \text{ TCID}_{50}/\text{ml}$
Human Coronavirus	229E	$1 \times 10^{1.5} \text{ TCID}_{50}/\text{ml}$
	OC43	$1 \times 10^4 \text{ TCID}_{50}/\text{ml}$
	NL63	$1 \times 10^4 \text{ TCID}_{50}/\text{ml}$


There was no cross-reaction and interference with the potential cross-reacting microorganisms listed below except SARS-CoV.

[Specimens in transport media]

1. Using a pipette, collect the 350μl of specimen from the collection cup or VTM. Mix the specimen with an extraction buffer.
2. Press the nozzle cap tightly onto the tube.


3. Few times mix







4. Rotate the swab 3-4 times against the nasopharyngeal wall.

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