

O.K.N.V.I . RESIST-5

Detect and identify the major Carbapenemases from 2 types of samples and Bacterial Culture

Flexible



✓ Positive Blood culture



✓ Rectal Swab

✓ Bacterial culture



No equipment needed !



ReSCape

NEW



For the screening of CP-BGN carriage by asymptomatic patient

- ✓ Reduced turn-around time (TAT)
- ✓ Better infection control management in outbreak setting

RESIST-BC


NEW




Fasten time of detection of CP-GNB in patients with bloodstream infection

- ✓ Faster medical decisions
- ✓ Help clinicians to optimize patient management


ReSCape



RESIST-BC



Save precious time, compared to conventional methods !!



Products	Targets	Description
Performance comparable co molecular methods		
OXA-48 K-SET	OXA-48-like	20 cassettes, buffer, 20 tube-droppers
KPC K-SET	KPC	20 cassettes, buffer, 20 tube-droppers
RESIST-3 O.O.K. K-SET	OXA-48-like, OXA-163*, KPC	20 cassettes, buffer, 20 tube-droppers
RESIST-3 O.K.N. K-SET	OXA-48-like, KPC, NDM	20 cassettes, buffer, 20 tube-droppers
OXA-23 K-SET	OXA-23-like	20 cassettes, buffer, 20 tube-droppers
O.K.N.V.I. RESIST-5	OXA-48-like, KPC, NDM, VIM, IMP	2x20 cassettes, buffer, 20 tube-droppers
RESIST-BC	OXA-48-like, KPC, NDM, VIM, IMP Not yet validated with OXA-163* ans OXA-23	Reagents kit for blood culture 1 RBCL solution, 1 MS buffer, 1 Washing buffer
ReSCape	OXA-48-like, OXA-163*, OXA-23, KPC, NDM, VIM, IMP	Culture reagent kit for CP-GNB on rectal swab 20 CProbeBe_broth tubes, 2 antibiotics vials, 2 sterile water tubes

Distributed by:

bestbion dx GmbH
Horbeller Str. 33
50858 Köln
Deutschland
Telefon: +49 2234 98795 - 0
Telefax: +49 2234 98795 - 29
Email: service@bestbion.com
Internet: www.bestbion.com

 **bestbion**^{dx}
smarter diagnostics

MADE IN BELGIUM

D-CUS-10-RES/01
October 2021